


INSTRUCTION FOR USE
VANTELIN COMPRESSION WEAR
Calf Sleeve

SIZE GUIDE

M: 33-38 cm

L: 38-43 cm

Measure the circumference (cm) around the thickest part of the calf as shown in the illustration.



FEATURES

Features

X-shaped taping line lifts calf muscle to suppresses excessive movement causing fatigue.

Highly breathable, all-day comfort

Highly absorbent material enhances wearing comfort by quickly absorbing and diffusing sweat.

Easy to stretch and move
with highly stretchable LYCRA® Fiber applicable for wide range of movement.

UV Protection

A 3D cutaway diagram of a calf sleeve. It shows an X-shaped taping line on the back of the calf. Red arrows point upwards from the taping line, indicating muscle lifting. Labels 'Upper anchor' and 'Taping line' are connected to their respective parts of the sleeve. The sleeve is shown over a foot and ankle.

INDICATION

Reduce excessive muscle movement causing fatigue.

INSTRUCTION FOR USE

1. Slide leg into calf sleeve with bottom of the "V" mark pointing down.
2. Align with the heel and pull up.
3. As shown in the illustration, adjust the position so that the calf muscle is supported from below and there is no slack. (Avoid covering knee joint.).



Wear directly on the skin. When wearing socks, wear them over the product. M size has a mark (II) and L size has a mark (III) in the same color as the sleeve, inside the upper anchor.

WARNING

- Do not use if you have: (1) abnormalities including swelling and strong pain due to sprain, fracture, or ligament injury etc., (2) peripheral vascular disorder, (3) abnormality such as injury, swelling, rash at site where support is to be worn, or (4) chronic dermatitis / allergy to chemical textile or rubber. If you experience symptoms such as rash / redness, itchiness, eruption, congestion, abnormal pain, or swelling, stop using immediately and consult your healthcare professional. Stop using the product immediately if you feel discomfort.
- Use only as directed.
- Use correct size.
- Do not use when sleeping.
- Do not wear two layers.
- It may cause blood circulation disorder, skin disorder, or neurological disorder, put on and take off the product repeatedly as appropriate when using for a long time.
- * If the product is not positioned correctly, it may not function as intended.
- * Functions obtained by wearing the product may vary depending on each person.

CARE INSTRUCTIONS (Use laundry net)

1. Wash separately below 40°C.
2. Do not bleach.
3. Do not tumble dry.
4. Hang and air-dry. Avoid direct sunlight.
5. Do not iron.
6. Do not dry clean.



PRECAUTIONS CONCERNING STORAGE AND HANDLING

1. Do not place the product close to fire and heat.
2. Keep out of reach of children.
3. There is a risk of the product tearing if worn in a wrong way or strong force is applied to part of the product.

PACKAGING

2 pieces / unit box

For both right and left calves

PRODUCT OWNER

Kowa Company, Ltd.

4-14, Nihonbashi-honcho 3-chome, Chuo-ku, Tokyo, JAPAN